

## ZOLL Defibrillator Competency Checklist

Please use this checklist when "playing" with the machine. I have left some reference sheets to review if You have questions, plus I will do roving live in-services as well.

<b>MONITOR</b>	
Locate pediatric paddles ➤ Able to remove and reattach paddles to the adult paddles	
Turns monitor ON	
Changes Leads	
Adjust size of waveform	
Locate Volume control key Turn QRS tone off ➤ QRS tone always defaults to ON	
Locate alarm button Change upper and lower limits of an alarm Suspend an audio alarm ➤ Alarm default is to OFF	
<b>DEFIBRILLATOR</b>	
Turns defibrillator ON	
Locates joule indicator ➤ Adjusts joules ➤ Pediatric dosing for defibrillation and Cardioversion is <b>UNCHANGED</b>	
Locates synchronized button	
Locates charge button	
Locates discharge button	
<b>PACING</b>	
Place in pacing mode	
Adjust ma and rate	
Locate multifunction pad, place on machine	
<b>RECORDING</b>	
Stop paper recording Push twice to stop paper recording	
Locate record key Push to receive paper recording of event and to mark a critical event	
<b>TESTING</b>	
Testing: Performed at 30J Test into the machine	
<b>OTHER</b>	
Change Paper: ➤ Slide out paper drawer ➤ Insert paper in the direction of black arrows ➤ The paper must touch rollers and close drawer ➤ Paper provided is square 2x2 size	
Summary key ➤ Locate print chart key ➤ (You can print all or a range) ➤ Locate print log (snapshot of all critical events EXCEPT ECG)	
<b>BATTERY</b>	
Remove a battery Locate battery indicator lights Locate recalibration indicator ➤ Batteries will trickle charge when plugged in ➤ How long the battery lasts depends on the type of use ➤ 35 defibrillations, 1.5 hours continuous ecg, 1.0hour of ecg and pacing	