

ZOLL Defibrillator Competency Checklist

Please use this checklist when "playing" with the machine. I have left some reference sheets to review if You have questions, plus I will do roving live in-services as well.

MONITOR	
Locate pediatric paddles ➤ Able to remove and reattach paddles to the adult paddles	
Turns monitor ON	
Changes Leads	
Adjust size of waveform	
Locate Volume control key Turn QRS tone off ➤ QRS tone always defaults to ON	
Locate alarm button Change upper and lower limits of an alarm Suspend an audio alarm ➤ Alarm default is to OFF	
DEFIBRILLATOR	
Turns defibrillator ON	
Locates joule indicator ➤ Adjusts joules ➤ Pediatric dosing for defibrillation and Cardioversion is UNCHANGED	
Locates synchronized button	
Locates charge button	
Locates discharge button	
PACING	
Place in pacing mode	
Adjust ma and rate	
Locate multifunction pad, place on machine	
RECORDING	
Stop paper recording Push twice to stop paper recording	
Locate record key Push to receive paper recording of event and to mark a critical event	
TESTING	
Testing: Performed at 30J Test into the machine	
OTHER	
Change Paper: ➤ Slide out paper drawer ➤ Insert paper in the direction of black arrows ➤ The paper must touch rollers and close drawer ➤ Paper provided is square 2x2 size	
Summary key ➤ Locate print chart key ➤ (You can print all or a range) ➤ Locate print log (snapshot of all critical events EXCEPT ECG)	
BATTERY	
Remove a battery Locate battery indicator lights Locate recalibration indicator ➤ Batteries will trickle charge when plugged in ➤ How long the battery lasts depends on the type of use ➤ 35 defibrillations, 1.5 hours continuous ecg, 1.0hour of ecg and pacing	